



Clorox Commercial Solutions®

**Clorox® Germicidal Wipes: 150 ct, 6"x5"**

Your Solution for F Tag 441 Disinfection Compliance



### Ensure F Tag 441 Compliance:

The Centers for Medicare & Medicaid services F Tag 441 document indicates that an Immediate Jeopardy condition exists any time in which a glucometer is not properly disinfected between uses from one individual to another, and also indicates that surveyors will be evaluating whether this is taking place. Additionally, some state organizations now require the use of bleach to disinfect these items.

**Clorox® Germicidal Wipes are your solution to three important issues that state surveyors will observe in your facility's cleaning and disinfecting practices:**

- ✓ Small non-disposable equipment is properly cleaned and disinfected after each use.
- ✓ High-touch surfaces are cleaned and disinfected.
- ✓ Equipment in transmission-based precaution rooms are thoroughly cleaned and disinfected.

### The perfect solution:

Incorporate Clorox® Germicidal Wipes into your facility's standard precautions and help protect those under your care. Clean and disinfect your healthcare equipment regularly with 150 count Clorox® Germicidal Wipes, now available in a convenient 6"x5" size for small surface disinfection.



### Clorox® Germicidal Wipes are:

- Designed for a fast-paced environment: kills many pathogens in one minute or less
- EPA and state registered to kill 36 different pathogens:
- Formulated with a 5500 ppm concentration of Clorox® bleach
- Engineered with odor-masking agents
- Tested for a 1-year shelf life
- Premoistened to ensure proper concentration: helps protect patients and staff from spills & chemical odors

Contact Time	Bacteria	<b>30 sec</b>
	Viruses (Including Norovirus)	<b>1 min</b>
	Tuberculosis	<b>2 min</b>

See product label for specific organisms

### Use with confidence:

Clorox® Germicidal Wipes contain anti-corrosion agents and have been tested for use on many common non-porous medical equipment surfaces, including plastics, glass, healthcare-grade aluminum, stainless steel, synthetic rubber, and more. Avoid contact with surfaces such as wood, painted surfaces, paper, and cloth.

### Where else to use Clorox® Germicidal Wipes:

In addition to glucometers, Clorox® Germicidal Wipes are designed for use on non-critical medical devices\* such as:

- Blood pressure cuffs
- Stethoscopes
- Ultrasound equipment
- Over-bed tables
- Other high-touch surfaces

For more information, visit [www.cloroxprofessional.com](http://www.cloroxprofessional.com), email [healthcare@clorox.com](mailto:healthcare@clorox.com), or contact your Clorox representative for a demonstration.

\*Always use product as directed on hard, non-porous surfaces.

If unsure about medical equipment surface compatibility, please follow equipment manufacturer guidelines.

Centers for Medicare & Medicaid Services State Operations Manual, Appendix PP, F Tag 441 p. 38

© 2009 Clorox Professional Products Company, 1221 Broadway, Oakland, CA 94612, [www.cloroxprofessional.com](http://www.cloroxprofessional.com)